


Installation Instructions Padded Seat Option for Rival™ Athletic Locker

Basic Instructions to be used with Rival Athletic Lockers only!
For more detailed product information, refer to the
Rival Athletic Locker Installation Instructions 243A242.

1. For better adhesion, clean the locker seat with the included alcohol wipe pads.
2. Peel the paper from the hook and loop strips at the bottom of the seat.
3. With the seam facing the rear, center the padded seat to be flush with the front of the locker seat and push down to adhere.



Visit gearboss.com/rival-lockers.php for detailed instructions and videos.

Note: Please read and understand these instructions before assembling.

Note: If you need additional information, contact Wenger Corporation using the information below.